

# Hints to familiarise yourself with the N6 sound processor before 'switch-on'



## Sound processor

### How do I turn my processor on?

Press and hold the **lower button** for three seconds, or simply attach a charged battery.

### How do I turn my processor off?

Press and hold the **lower button** for three seconds, or remove the battery.

### Why do I have different programs?

Often in the first month or so, as you get used to the sound, the clinician may give you up to four programs of different volumes or settings.



## How do I change programs?

This can be done via either the sound processor or the remote assistant:

### (a) Via the sound processor:

Quick press of **lower button** for one second. The indicator light will flash green, corresponding to the program number you are on. For example two green flashes indicate you are on program number two.

### (b) Via the remote assistant:

Press the 'Change Program' button that corresponds to the program number/listening environment you would like. The remote assistant can also be used for trouble-shooting.



## Batteries

### Changing rechargeable batteries?

Twist battery module to the left to separate from the sound processor. Twist a charged battery module to the right to attach to the processor.



### How do I change the disposable batteries?

Slide cover down to replace the batteries. Insert batteries positive (+) side up. Keep the battery module attached to processing unit when changing batteries.



## Who to contact if you have a concern

For any concerns, please contact:

### Cochlear Implant Clinic

6th Floor, Smorgon Family Wing

32 Gisborne Street

East Melbourne VIC 3002

Ph: 03 9929 8624

For medical concerns, please visit the Eye and Ear Emergency Department or call 03 9929 8666.

**Disclaimer** This document describes the generally accepted practice at the time of publication only. It is only a summary of clinical knowledge regarding this area. The Royal Victorian Eye and Ear Hospital makes no warranty, express or implied, that the information contained in this document is comprehensive. They accept no responsibility for any consequence arising from inappropriate application of this information.

Hints to familiarise yourself with the sound process before switch on #203 | Owner: Cochlear Implant Clinic | Last published: 17/02/2016 | Next review: 17/02/2018

